

Pragma: no-cache  
X-Content-Type-Options: nosniff  
X-XSS-Protection: 1; mode=block  
X-Frame-Options: DENY  
Ocp-Apim-Trace-Location:  
https://apimgmtstln1vrkzf9gvrl3n.blob.core.windows.net/apiinspectorcontai  
ner/gthPseliu-KjMEue3RdyCg2-102?sv=2015-07-  
08&sr=b&sig=CYgXrJXEPz4l3aNWn1WbH9VQGjnWivT1N5dvu0uXACA%3D&se=2017-06-  
22T11%3A21%3A15Z&sp=r&traceId=03fe0c468ae24bleb75de6cdda481f94  
Cache-Control: no-store, must-revalidate, no-cache, max-age=0  
Date: Wed, 21 Jun 2017 11:21:15 GMT  
Content-Length: 8247  
Content-Type: application/json; charset=UTF-8  
Expires: 0

```
[
  {
    "ChemistryReportEHRExtract": [
      {
        "PerformingUnit": {
          "UnitID": {
            "extension": "CC000123",
            "root": "1.2.752.129.2.1.4.1"
          },
          "UnitLocalID": {
            "extension": "",
            "root": "LocalUnitOID"
          }
        },
        "ehr_id": {
          "extension": "SE162321000024-  
Cambio_Common_U:EhrId:198611062384",
          "root": "1.2.752.129.2.1.2.1"
        },
        "Report": {
          "ReportDateTime": {
            "value": "20170621132032"
          },
          "ReportId": "",
          "TestResult": [
            {
              "TestId": "NPU02593",
              "TestName": "B-Leukocyter",
              "Pathological": "0",
              "IUPAC": "",
              "Unit": "*10^9 /L",
              "Range": "3,5 - 8,8 *10^9 /L",
              "rc_id": {
                "extension": "SE162321000024-  
Cambio_Common_U:RcId:ChemistryResult:84346",
                "root": "1.2.752.129.2.1.2.1"
              },
              "Result": "6,2"
            },
            {
              "TestId": "NPU01807",
              "TestName": "P-Kreatinin",
              "Pathological": "0",
            }
          ]
        }
      }
    ]
  }
]
```

```

        "IUPAC": "",
        "Unit": "µmol/L.",
        "Range": "45-90",
        "rc_id": {
            "extension": "SE162321000024-
Cambio_Common_U:RcId:ChemistryResult:84343",
            "root": "1.2.752.129.2.1.2.1"
        },
        "Result": "83"
    }
],
"ReportStatus": "CO",
"rc_id": {
    "extension": "SE162321000024-
Cambio_Common_U:RcId:ChemistryReport:98048",
    "root": "1.2.752.129.2.1.2.1"
},
"ReportVersionDateTime": {
    "value": "20170621132032"
}
},
"ehr_system": {
    "extension": "SE162321000024-Cambio_Common_U",
    "root": "1.2.752.129.2.1.2.1"
},
"Request": {
    "RequestId": "80000461",
    "CollectionDateTime": "20170521140100",
    "RequestDateTime": {
        "value": "20170621131700"
    },
    "Priority": "Rutin",
    "rc_id": {
        "extension": "SE162321000024-
Cambio_Common_U:RcId:ChemistryRequest:84329",
        "root": "1.2.752.129.2.1.2.1"
    },
    "RequestCommentary": "Nedsatt hörsel hos patienten. Patienten
är också i behov av rullstol och betalar i Pund (£) . Provsvaret kommer
ha enheten µmol/L."
},
"ReceivingUnit": {
    "UnitID": {
        "extension": "SE5565189692-101JW",
        "root": "1.2.752.129.2.1.4.1"
    },
    "UnitLocalID": {
        "extension": "",
        "root": "LocalUnitOID"
    }
},
"Patient": {
    "PatientID": {
        "extension": "198611062384",
        "root": "1.2.752.129.2.1.3.1"
    },
    "PatientProtectedId": false,
    "PatientName": {

```

```

    "PatientFamilyName": "Alm",
    "PatientMiddleName": "",
    "PatientGivenName": [
      {
        "qualifier": "NN",
        "value": "Ulla"
      }
    ]
  },
  "PatientInternalID": "4252"
},
"CareProviders": {
  "CareProvider": [
    {
      "CareProviderUnit": {
        "UnitID": {
          "extension": "Not implemented",
          "root": "1.2.752.129.2.1.4.1"
        },
        "UnitLocalID": {
          "extension": "Not implemented",
          "root": "1.2.752.129.2.1.4.1"
        }
      },
      "CareProviderInternalID": "39747",
      "ProviderPosition": {
        "CareProviderPositionID": {
          "extension": "Not implemented",
          "root": "Define"
        },
        "CareProviderPositionName": "Läkare"
      },
      "CareProviderType": "RequestingPhysician",
      "CareProviderID": {
        "extension": "",
        "root": "1.2.752.129.2.1.4.1"
      },
      "CareProviderName": "Sven Svensson"
    },
    {
      "CareProviderUnit": {
        "UnitID": {
          "extension": "Not implemented",
          "root": "1.2.752.129.2.1.4.1"
        },
        "UnitLocalID": {
          "extension": "Not implemented",
          "root": "1.2.752.129.2.1.4.1"
        }
      },
      "CareProviderInternalID": "39747",
      "ProviderPosition": {
        "CareProviderPositionID": {
          "extension": "Not implemented",
          "root": "Define"
        },
        "CareProviderPositionName": "Läkare"
      }
    }
  ]
}

```

```

    "CareProviderType": "ReportReceivingPhysician",
    "CareProviderID": {
      "extension": "",
      "root": "1.2.752.129.2.1.4.1"
    },
    "CareProviderName": "Sven Svensson"
  }
]
},
"Units": {
  "Unit": [
    {
      "UnitID": {
        "extension": "SE162321000024-0001",
        "root": "1.2.752.129.2.1.4.1"
      },
      "UnitName": "Landstinget i Länet",
      "UnitLocalID": {
        "extension": "",
        "root": "LocalUnitOID"
      }
    }
  ],
  {
    "PhoneNumbers": {
      "Phone": [
        {
          "PhoneComment": "",
          "PhoneType": "Växel",
          "PhoneNumber": "023-490000"
        }
      ]
    }
  },
  "UnitAddresses": {
    "UnitAddress": [
      {
        "UnitAddressCountry": "Sverige",
        "UnitAddressPostcode": "79182",
        "UnitAddressAddress": "KIR mott- Falu Lasarett",
        "UnitAddressCity": "Falun",
        "UnitAddressCounty": "Dalarna",
        "UnitAddressComment": "",
        "UnitAddressType": "Besöksadress"
      }
    ]
  },
  "UnitID": {
    "extension": "SE5565189692-101JW",
    "root": "1.2.752.129.2.1.4.1"
  },
  "UnitName": "Kirurgkliniken",
  "EmailAddresses": {
    "Email": [
      {
        "EmailType": "",
        "EmailAddress": "info@kirmott.se"
      }
    ]
  }
},

```

```

"UnitSecrecyGroups": {
  "Group": [
    "Kirurgkliniken"
  ]
},
"HeadCareproviderUnitID": {
  "extension": "SE162321000024-0001",
  "root": "1.2.752.129.2.1.4.1"
},
"UnitLocalID": {
  "extension": "",
  "root": "LocalUnitOID"
}
},
{
  "PhoneNumbers": {
    "Phone": [
      {
        "PhoneComment": "",
        "PhoneType": "HSA_Växeltелефон",
        "PhoneNumber": "0123-456789"
      }
    ]
  },
  "UnitAddresses": {
    "UnitAddress": [
      {
        "UnitAddressCountry": "Sverige",
        "UnitAddressPostcode": "98700",
        "UnitAddressAddress": "Labbvägen 23",
        "UnitAddressCity": "emstad",
        "UnitAddressCounty": "Labblän",
        "UnitAddressComment": "",
        "UnitAddressType": "HSA_Besöksadress"
      }
    ]
  },
  "UnitID": {
    "extension": "CC000123",
    "root": "1.2.752.129.2.1.4.1"
  },
  "UnitName": "Kemlabb AB",
  "EmailAddresses": {
    "Email": [
      {
        "EmailType": "HSA_Email",
        "EmailAddress": "info@kemlabb.se"
      }
    ]
  },
  "UnitLocalID": {
    "extension": "",
    "root": "LocalUnitOID"
  }
}
]
},
"HeadCareproviderUnit": {

```

```
"UnitID": {
  "extension": "SE162321000024-0001",
  "root": "1.2.752.129.2.1.4.1"
},
"UnitLocalID": {
  "extension": "",
  "root": "LocalUnitOID"
}
},
"HCMContact": {
  "HCMContacInternaltID": "97983",
  "target_id": {
    "extension": "SE162321000024-
Cambio_Common_U:RcId:Contact:97983",
    "root": "1.2.752.129.2.1.2.1"
  },
  "HCMContactBegin": {
    "value": "20170517230708"
  },
  "HCMContactType": "Mottagningsbesök",
  "ContactSecrecyLevels": [
    {
      "SecrecyLevel": ""
    }
  ],
  "HCMContactStatus": "Ankommen",
  "ContactPerformingUnit": "1062"
},
"RequestingUnit": {
  "UnitID": {
    "extension": "SE5565189692-101JW",
    "root": "1.2.752.129.2.1.4.1"
  },
  "UnitLocalID": {
    "extension": "",
    "root": "LocalUnitOID"
  }
},
"time_created": {
  "value": "20170621132040"
},
"AccessZoneUnit": {
  "UnitID": {
    "extension": "",
    "root": "1.2.752.129.2.1.4.1"
  },
  "UnitLocalID": {
    "extension": "",
    "root": "LocalUnitOID"
  }
}
}
]
}
```